



Collection of Bag Urine from a Baby

Please follow these instructions:

1. Bath baby washing the penis/labia well.
2. Pat the baby dry. Do not apply cream or powder.
3. Remove the paper backing to expose the adhesive.




4. Place the opening of the bag so that it covers the upper half of the genitals (Diagram 2) for a girl, or, for a boy place the penis through the opening (Diagram 3).




5. Look at the bag approximately half hourly until there is urine present.

6. Remove the bag from the baby and fold the sticky surfaces together.



7. Open the yellow top specimen jar being careful not to touch the inside of the jar or lid. Clean the bottom of the bag with the alcohol wipe. Wait for 2-3 minutes then using scissors snip a bottom corner so that the urine drains into the jar.



8. Screw the lid firmly onto the jar and write baby's name, date of birth, date / time of collection on the label.
9. Put container into your fridge until you can deliver to your nearest Western Diagnostic Pathology Collection Centre.
10. If you have any queries, please telephone us on **9317 0999** and speak to the Duty Manager.